




**Solutions for Home Sharps Collections**

CPSF Forum – January 18, 2012



**Disposal Solutions for Home Users**




Mailback Kits      Community Collection




**Mailback Kits**

A **SAFE** and **CONVENIENT SOLUTION** for the disposal of home sharps.

Easy system includes all components required for collection and transportation of needles, lancets, or syringes.

- Step-by-step instructions
- Puncture-resistant sharps container with leak proof screw cap
- Bar-coded shipping containers
- Prepaid and pre-addressed shipping labels
- Tracking documents to audit the disposal process

1.4 Quart (101V4)  
3.34" x 3.102" x 8"  
Holds approximately  
45 3cc syringes



**Mailback Kits**

- System includes all components for collection and transportation of sharps via US Postal Service and disposal by Stericycle
- Pre-paid and pre-addressed shipping labels add ease and convenience for end users
- Tracking documents audit shipping and disposal process and assure the generator of proper disposal
- Available in bulk quantities for a central distribution point in the community or can be set-up for direct shipment to the resident



**Community Collection**

The **safe, affordable** solution for community sharps disposal





**Community Collection**

- Designed specifically for hospitals, pharmacies, and local government buildings
- Collection location purchases the kiosk, replacement internal sharps containers and waste removal services
- Collection location's staff routinely manages the waste, by replacing full containers with new, empty containers
- Requires a regularly scheduled pick-up service



**Stericycle**  
Protecting People. Reducing Risk.

## Let's Talk About Design...



**DURABLE**  
Heavy-gauge steel cabinet holds up to daily public use.

**SAFE**  
Secure locking mechanism prevents tampering.

Dimensions:  
43" tall x 22" wide x 18" deep

Easy-to-read instructions provided in English and Spanish

Special drawer design prevents access to interior container

Easy-glide drawer safely accepts sharps containers as large as 7 gallons

7-gallon sharps container collects sharps and needles—fits securely for secure sharps disposal

7

**Stericycle**  
Protecting People. Reducing Risk.

## Instructions for Users

**WARNING**  
PERILOSO

**BIOHAZARD**

**Stericycle**  
Needle Drop Box

### HOME SHARPS DISPOSAL ONLY

PARA LA DISPOSICION DE CASA AGUDO SOLO

- 1 Pull handle to fully open drop box drawer.**  
Jale la manijeta completamente para abrir la gaveta de depósito.
- 2 Place sealed sharps container inside.**  
Coloque adentro, completamente cerrado, los recipientes de "sharps".
- 3 Slowly close drop box drawer.**  
Cierre la caja lentamente.

**DO NOT** dispose of pharmaceutical or hazardous waste (old medicines, batteries, motor oil, paint, etc.)  
NO - depositar desperdicios peligrosos o farmacéuticos (medicinas viejas, baterías, aceite de auto, pintura, etc.)

8

**Stericycle**  
Protecting People. Reducing Risk.

## Thank You

Speaker Contact:  
John Katchka  
630-797-7300  
jkatchka@stericycle.com

9